




Your local partner for all your orders
Chemical products - Paper or plastic consumables - Glassware - Apparatus



A market leader in labeling and
detection systems for biological
science



Comprehensive offering across the IHC and IF workflow
Very responsive and well-documented support on their website

Reactolab offers 30% discount on their catalog

BONUS: the latest from Vector !

VectaMOUNT EXPRESS Mounting Medium

- Non-aqueous cleaning and mounting medium for IHC-stained slides
- Mount after 2 min. washing !
- Save dehydrating and cleaning solvents.
- Optimal viscosity and refractive index (1.49)
- 60ml : app. 1200 tissue sections

VECTASHIELD PLUS Antifade Mounting Media

- New formulation with reliable signal retention of all commonly used fluorophores
- Sections can be viewed in 1 hour after mounting
- With or without DAPI counterstain
- Ready to use

H-5700-60: 60ml: ~~143~~ 80 frs HT

H-1900-2 (-DAPI) : 2ml: ~~81~~ 45 frs HT
H-2000-2 (+DAPI): 2ml: ~~85~~ 45 frs HT